

PNEUMONIA DETECTION

Deanna Gould

Data Science Flex

THE PURPOSE

- Predicting pneumonia with X-Ray images for HealthWorx, a telehealth company
- Decrease emergency room wait times
- Refer a patient to a third-party radiology company

```
train_set:
```

```
-----
```

```
Pneumonia = 3476
```

```
Normal = 942
```

```
val_set:
```

```
-----
```

```
Pneumonia = 409
```

```
Normal = 409
```



```
test_set:
```

```
-----
```

```
Pneumonia = 390
```


```
Normal = 234
```

PNEUMONIA OR NOT?




PIXEL COMPARISON

Pixel Intensity Distribution for a Normal X-Ray




Pixel Intensity Distribution for a Pneumonia X-Ray



FINAL CONFUSION MATRIX

- - 98.3% true positive rate
- - 1.7% false negative rate
- - 5.6% false positive rate



RECOMMENDATIONS

- - Use Inception V3 for modeling purposes
- - Generate patient data to ensure we're training the model on diverse demographics
 - - Age, sex, race, smoking history
- With more time, try other models, and with more resources, use more X-Ray images

QUESTIONS

If there are questions after the conclusion of the presentation, I can be reached at deannagould4@gmail.com.